



PAR-H50



Portarotoli

Rollholder

gandus®
SALDATRICI
PACKAGING MACHINES



Caratteristiche costruttive - *Technical features*

Portarotoli per buste di sterilizzazione con tagliatore bidirezionale, il **PR-H50** può essere appeso alla parete o appoggiato sopra un piano di lavoro in prossimità di una termosaldatrice per la realizzazione di buste, nella lunghezza desiderata, senza ricorrere a forbici o ad altri utensili. Costruzione in acciaio Inox AISI 304 e alluminio anodizzato.

*Roll dispenser for medical pouches, knife cuts in both directions, **PR-H50** can be mounted on a wall or placed on a working table close to the sealing machine. Dispenser is useful where different lengths of pouch are required without the need of scissors or any other cutting tool. Made in stainless steel AISI 304 and anodized aluminum.*

Dati Tecnici

Larghezza max di taglio	510 mm
Cursore bidirezionale	
Lama autoprotetta a 4 taglienti	
Dimensioni	630 x 280 x 360 mm
Peso	8,5 Kg

Technical Data

<i>Max. Cut length</i>	<i>510 mm</i>
<i>Bidirection knife</i>	
<i>Autosafe 4 cutting side-blade</i>	
<i>Dimensions</i>	<i>630 x 280 x 360 mm</i>
<i>Weight</i>	<i>8,5 Kg</i>

COMPANY
WITH QUALITY MANAGEMENT
SYSTEM CERTIFIED BY DNV
=ISO 9001:2000=

gandus[®]
SALDATRICI
PACKAGING MACHINES

